



LERNAPHARM

Specialized in the Development and the Production of Antiseptics and Disinfectants

● **Over A Quarter Century Of Service**

● **International Presence**

**PHARMACEUTICAL
MANUFACTURING**



**Promoting
Health And Wellness**





LORIS™ Medium Alcohol Swabs

Product Code : 103-01



Characteristics

How To Use

Tear open package and apply a wipe directly to the skin. Rub skin briskly in a circular motion from injection site outward.



2323 Halpern
St-Laurent (Montreal)
Quebec, Canada, H4S 1S3

Tel: + 514-331-4634
Toll-Free: 1-866-322-4634
Fax: + 514-331-5488
info@lernapharm.com
www.lernapharm.com








LORIS™ Medium Alcohol Swabs

Product Code : 103-01

Minimal Acceptance Criteria for Antiseptics Microbial Limit Test

Test	USP Nonsterile Dosage Forms(Cutaneous Use)	PVP-I Solution
Total Aerobic Microbial Count (TAMC)	≤100 cfu / mL	<10 cfu / mL
Total Combined Yeast and Molds Count (TYMC)	≤10 cfu / mL	<2 cfu / mL
taphylococcus aureus detection	Absent	Absent
Pseudomonas aeruginosa detection	Absent	Absent
B. Cepacia	Not specified	Absent

Packaging Information

Product	Quantity	Size(L,W,H)	Weight
Wipe 	-	120 mm x 220 mm	
Pouch 	1 x Wipe per Pouch	65 mm x 35 mm	
Display Box 	100 x Pouch per Display Box	11 cm x 7 cm x 6 cm	
Case 	40 x Display Box per Case	31 cm x 23 cm x 30 cm	5 kg
Pallet 	102 x Case per Pallet	122 cm x 102 cm x 192 cm	473 kg