



Technical Data Sheet

Product: Loris Alcohol Swabs

NPN 02248362 (Canada)

Product Packaging

Product Code: 103-03

Pouch size: 50 mm x 50 mm

Open wipe: 65 mm x 30 mm

Closed wipe: 32 mm x 30 mm

Seals dimensions: 3 mm all around

Solution content: Minimum 0.35 mL/pouch

Packaging Format: 200 pads / display box & 20 display boxes / corrugate

Corrugates per skid: 144

Corrugates per row: 18


Rows per skid: 8

Product Label

Use: Antiseptic skin cleanser prior to injection.

Directions: Tear open package and rub skin briskly in a circular motion from injection site outward.

Cautions and warnings: For external use only. Avoid contact with eyes. If contact occurs, flush eyes with water. Discontinue use and consult a health care practitioner if irritation develops. Store in a cool dry place and away from sources of ignition.





Technical Data Sheet

In-Process Quality Controls on Finished Products

Test	Method	Specification
Fill volume	S-MAN-003	Minimum 0.35 mL
Seal Integrity	S-MAN-013	No leaks found in bundle of pouches when hardly squeezed between two plates

Chemical / Biological Analysis

Test*	Method	Specification
Appearance	TL-S-QUA-010	Transparent, clear, colorless and clean liquid with a characteristic odor
Specific Gravity	TL-S-QUA-008	0.872-0.883 at 20°C
Assay of Isopropyl Alcohol by specific gravity determination	TL-S-QUA-023	68.0-72.0%
Total Aerobic Microbial Count (TAMC)	TL-S-QUA-028	< 10 cfu/ml
Total Combined Yeast and Mold Count (TYMC)	TL-S-QUA-028	< 2 cfu/ml
S. aureus	TL-S-QUA-053	Absence
P. aeruginosa	TL-S-QUA-053	Absence
B. cepacia	TL-S-QUA-053	Absence

*Tests are made on the bulk solution